

IMPROVING TECHNOLOGIES
IMPROVED HEALTH

# PRODUCT CATALOGUE

IMPROVING TECHNOLOGIES
IMPROVED HEALTH



EARLY JUDGEMENT of RELATED DISEASES RAPID INFECTION CLASSIFICATION ACCURATE TREATMENT EVALUATION



# COMPANY PROFILE



Established in Beijing, 2013, Pro-med has been qualified as National High-tech Enterprise based on professional IVD detecting service, which contains instruments and agents of POCT, coagulation and molecular diagnosis.

With deep concern about public healthcare and clinical demands, Pro-med develops different PMDT qualitative/quantitative analysis systems for diseases diagnosis such as respiratory pathogens detecting, bacteria/virus infection distinguishing, AMI & HF managing, gastric system disorder, low health level and series of common infectious diseases rapid testing.

As for clinical coagulation analysis, Pro-med provides an easy but clear solution by CWPS thrombelastograph, which could evaluate coagulation capacity through progress of thrombin activation, fibrin formation, platelet aggregation and fibrinolysis. Besides, a real-time coagulation graph could also be auto-drawn to visualize these progresses in the final report.

May the Pro-med products deliver pleasant experiences by accurate diagnosis, convenient operation and professional services to all users and clients!

### **IMPROVING TECHNOLOGIES, IMPROVED HEALTH**

# TABLE OF CONTENTS

01

# ACCURATE QUANTITATIVE ANALYSIS SYSTEM

PMDT9000 Immunofluorescence Analyzer—02 Feature Detection Kits—03



02

# RAPID QUALITATIVE ANALYSIS SYSTEM

Feature Colloidal Gold Detection Kits—08

03

### **THROMBELASTOGRAPH**

CWPS TEG analyzer—10
Feature Thrombelastograph Detection Kits—11



04

#### **MOLECULAR DIAGNOSTICS**


Real-time PCR system—10
Feature Molecular Detection Kits—11



# **01**ACCURATE QUANTITATIVE ANALYSIS SYSTEM

PMDT9000 Immunofluorescence Analyzer—02
Feature Detection Kits—03

# PMDT9000 IMMUNOFLUORESCENCE ANALYZER

















Power	AC220, 50Hz	
Detection time	5-15min	
Available value	0 ~ 70000mv	
Size	267*220*140mm	

# Feature Immunofluorescence Detection Kits

### **Thyroid Hormone**

·TSH Rapid Test

Thyroid function evaluation

·T3 Rapid Test

Specific diagnostic indicators of T3-type hyperthyroidism

·T4 Rapid Test

Sepcific diagnostic indicators of T4-type hyperthyroidism

+F13

Specific diagnostic indicators of FT3-type hyperthyroidism

·FT4

Specific diagnostic indicators of FT4-type hyperthyroidism



#### **Hormone**

· **β-HCG Rapid Test**Assist in judging pregnancy

· AMH Rapid Test

Ovarian health evaluation

· INH-B Rapid Test

Evaluation of testicular spermatogenesis function



### **Cardiac & Cerebral Diseases**

·ST2/NT-proBNP Rapid Test

Comprehensive management of heart failure

·ST2 Rapid Test

Early recognition; risk evaluation Prognostic assessment of heart failure

D-Dimer Rapid Test

Clotting condition determination

·NT-proBNP Rapid Test

Diagnosis of heart failure; Identification of dyspnea

· cTnI/Myo/CK-MB Rapid Test

Diagnosis of myocardial injury Evaluation of myocardial infarction ·S100 Rapid Test

Markers of brain injury

· CK-MB Rapid Test

Diagnosis of myocardial injury Evaluation of myocardial infarction

· B-type natriuretic Rapid Test

Diagnosis of heart failure Identification of dyspnea

·cTnl Rapid Test

Markers of myocardial injury

· Myo Rapid Test

Diagnosis of myocardial injury Evaluation of myocardial infarction



#### **Inflammation & Infection**

· Procalcitonin Rapid Test (PCT)

Systemic bacterial infections; Sepsis recognition

·C Reactive Protein (hs-CRP+CRP) Rapid Test

Assists in the diagnosis of sepsis and severe bacterial infections; Guidance on the rational use of antibiotics Assists in the differential diagnosis of local and systemic bacterial infections

· Serum Amyloid A Protein/C Reactive Protein (hs-CRP+CRP) Rapid Test (SAA/CRP)

Assists in identifying bacterial/viral infections; Auxiliary assessment of the course of infectious disease Guidance on the rational use of antibiotics


· Serum Amyloid A Protein Rapid Test (SAA)

Assists in identifying bacterial/viral infections; Auxiliary assessment of the course of infectious disease Guidance on the rational use of antibiotics

· IL-6 Rapid Test

Quick inflammatory determination

# Feature Immunofluorescence Detection Kits



### **Gastric Health**

- · Pepsinogen I, Pepsinogen II Rapid Test PG I /PG II Gastric mucosal injury judgment
- Gastrin-17 Rapid Test G17 Suggests gastric acid secretion; Gastric health testing indicators

### **Tumor**

- · CEA Rapid Test
  Assist in the diagnosis of malignancy
- · AFP Rapid Test
  An important marker of liver cancer prognosis



### Healthcare

- ·25-OH-VD
- Assess VD levels; Bone, immunity, cell growth related
- ·VB12 Rapid Test
- Nutrition level evaluation
- · N-MID Rapid Test
- Osteoporosis risk screening; Osteoporosis type discrimination/Endocrine-related
- · Ferritin Rapid Test
- Predicts iron deficiency anemia

## **PMDT 9000 CATALOG**

SERIES	PRODUCT NAME	REF	SAMPLE	CUT OFF	LINEAR RANGE
_, .,	TSH	PMDIA	serum/plasma/whole blood	0.3-5.6 μIU/mL	0.10~100.00μIU/mL
	T3	PMDIB	serum/plasma/whole blood	1.23~3.07 nmol/L	0.61~9.22 nmol/L
Thyroid Hormone	T4	PMDIC	serum/plasma/whole blood	66~181 nmol/L	12.87~300 nmol/L
Hormone	FT3	PMDID	serum/plasma/whole blood	3.1~6.8 pmol/L	0.40~50 pmol/L
	FT4	PMDIE	serum/plasma/whole blood	12~22 pmol/L	0.30~100 pmol/L
	β-HCG	PMDJA	serum/plasma/whole blood	various by age	5.00~50000.00 IU/L
Hormone	AMH	PMDJI	serum/plasma/whole blood	various by gestational week	0.1~18.00ng/mL
	INH-B	PMDJL	serum/plasma/whole blood	Male: 16.61~278.86pg/mL Female: various by age	10.00~900.00pg/mL
	sST2/NT-proBNP	PMDFE	serum/plasma/whole blood	sST2:<35.00ng/mL NT-proBNP:<300.00 pg/mL	sST2: 5.00~200.00ng/mL NT-proBNP: 100.00~20000.00pg/mL
	sST2	PMDFD	serum/plasma/whole blood	<35.00ng/ mL	5.00~200.00ng/mL
	NT-proBNP	PMDFB	serum/plasma/whole blood	<300.00 pg/mL	100.00~20000.00pg/mL
Cardiac	BNP	PMDFC	serum/plasma/whole blood	<100.00 pg/ml	50.00~3000.00pg/mL
&	cTnI/CK-MB/Myo	PMDFI	serum/plasma/whole blood	cTnl:<1.00ng/mL Myo:<70.00ng/mL CK-MB:<5.00ng/mL	cTnl: 0.30~25.00ng/mL Myo: 25.00~500.00ng/mL CK-MB: 2.50~300.00ng/mL
Cerebral	D-Dimer	PMDFA	serum/plasma/whole blood	<500ng/mL	200.00~5000.00ng/mL
Diseases	S100-β	PMDFH	serum/plasma/whole blood	<0.11ng/mL	0.08~10.00ng/mL
	cTnI	PMDFJ	serum/plasma/whole blood	<1ng/mL	0.2~40.00ng/mL
	CK-MB	PMDFL	serum/plasma/whole blood	<5.0 ng/ mL	3.00~80.00ng/mL
	MYO	PMDFK	serum/plasma/whole blood	<70.00ng/mL	50.00~400.00ng/mL
Inflammation	SAA/CRP	PMDGC	serum/plasma/whole blood	CRP:<10.00 mg/L hs-CRP:<3.00 mg/L SAA:<10.00 mg/L	SAA: 5.00~150.00mg/L CRP: 0.50~200.00mg/L
	SAA	PMDGB	serum/plasma/whole blood	<10.00 mg/L	5.00~150.00mg/L
& Infection	CRP	PMDGA	serum/plasma/whole blood	CRP:<10.00 μg/mL hs-CRP:<3.00 μg/mL	0.50~200.00mg/L
meetion	PCT	PMDGD	serum/plasma/whole blood	<0.50 ng/mL	0.10~40.00ng/mL
	IL-6	PMDGE	serum/plasma/whole blood	<10pg/mL	3.00~4000pg/mL
Gastric Health	PGI/II	PMDHB	serum/plasma/whole blood	PGI:>50.00ng/ml PGII:<15ng/ml	PGI: 5.00~100.00ng/mL PGII: 1.00~50.00ng/mL
	G17	PMDHC	serum/plasma/whole blood	1.00-7.00pmol/L	1.00~20.00pmol/L
Tumor	AFP	PMDOA	serum/plasma	<10.0ng/mL	3.0~2000.0ng/mL
Tullion	CEA	PMDOB	serum/plasma	<5 ng/mL	Correlation Coefficient (R) is $\geq$ 0.99
	25-OH-VD	PMDKB	serum/plasma/whole blood	27.99~81.81ng/mL	5.00~100.00ng/mL
Healthcare	FER	PMDKA	serum/plasma/whole blood	>20.00ng/mL	10.00~100.00ng/mL
Heattiicale	N-MID	PMDKD	serum/plasma	4.00ng/ml-30.00 ng/mL	3.00~100.00ng/mL
	VB12	PMDKC	serum/plasma	200.00~800.00pg/mL	150.00~2000.00pg/mL







PMDT9100 High efficiency:200T/H

# **Q2**RAPID QUALITATIVE ANALYSIS SYSTEM

Feature Colloidal Gold Detection Kits/07

# Feature Colloidal Gold Detection Kits

## COVID-19



- · COVID-19 Antigen Rapid Detection Kit
- · COVID-19 IgM & IgG Antibody Detection Kit
- ·COVID-19 Neutralizing Antibody Rapid Detection Kit
- ·Influenza A+B and COVID-19 Antigen Combo Rapid Detection Kit

# **Respiratory Tract**

- ·Adenovirus IgM Antibody Rapid Test
- · Respiratory Syncytial Virus IgM Antibody Rapid Test
- ·Influenza A Virus IgM Antibody, Influenza B Virus IgM Antibody and Parainfluenza Virus IgM Antibody Combined Rapid Test
- ·Influenza A virus and Influenza B Virus Antigen Combined Rapid Test
- · Mycoplasma Pneumoniae IgG antibody and IgM Antibody Rapid Test
- · Chlamydia Pneumoniae IgG Antibody and IgM Antibody Rapid Test
- · Mycobacterium Tuberculosis Antibody Rapid Test Kit
- · Clostridium Difficile GDH+ToxinA+ToxinB Antigen Rapid Test Kit (Colloidal Gold)
- ·Influenza A, Influenza B, COVID-19, Respiratory Syncytial Virus and Adenovirus Antigen Combo Rapid Detection Kit (Colloidal Gold)
- · Rotavirus /Adenovirus Antigen Rapid Test (Immunochromatography)

## **Gastric Health**

- · Helicobacter Pylori IgG Antibody Rapid Test
- · HP Antigen Rapid Test

# **Monkeypox Virus**

 $\cdot \, \mathsf{Monkeypox} \, \mathsf{Virus} \, \mathsf{IgG/IgM} \, \mathsf{Rapid} \, \mathsf{Test} \, \mathsf{Kit}$ 

# 03 THROMBELASTOGRAPH

CWPS TEG analyzer/10 Feature Thrombelastograph Detection Kits/11

# **CWPS-8800 TEG Analyzer**

# **Integrated TEG analyzer**

- $\cdot \mathsf{Semi\text{-}automated}\,\mathsf{TEG}\,\mathsf{analyzer}$
- · Integrated analyzer and computer, PC free
- $\cdot {\it Capable for all TEG tests in a single 4-channeled \ analyzer}$



Power	AC220, 50Hz	
Detection time	≥30min	
channel number	4	
Size	348*355*525mm	

## Feature Thrombelastograph Detection Kits



### **CK Test**

### Thrombelastograph (General Cup)

- · general detection of coagulation capability
- · risk evaluation of thrombus diseases and bleeding
- · stable QC for products quality

SAMPLE: CITRATED WHOLE BLOOD SPECIFICATION: 1/5/10/20/25/30/40/50 TEST/KIT







### AA, ADP, AA&ADP Test

#### **Platelet Aggregation Function Test**

- $\cdot \, \mathsf{Effect} \, \, \mathsf{evaluation} \, \, \mathsf{of} \, \, \mathsf{platelet} \, \mathsf{aggregation} \, \,$
- · Function test of anti-platelet drugs
- · AA for aspirin, ADP for clopidogrel (Plavix), tegrello and so on

SAMPLE: CITRATED&HEPARINED WHOLE BLOOD SPECIFICATION: 1/5/10/20/25/30/40/50 TEST/KIT



#### **CRT Test**

### **Activation of Blood Coagulation Test**

- · Rapid detection of coagulation capability
- · Fast coagulation judgement for emergency
- · Monitoring in surgery, ICU or even field trauma

SAMPLE: CITRATED WHOLE BLOOD SPECIFICATION: 1/5/10/20/25/30/40/50 TEST/KIT



### **CKH Test**

#### Thrombelastograph (Heparinase Cup)

- $\cdot \mathsf{Fast}\ \mathsf{detection}\ \mathsf{of}\ \mathsf{heparin}\ \mathsf{in}\ \mathsf{blood}$
- · Effect evaluation of heparin neutralization therapy
- ·Test to find real reason of bleeding

SAMPLE: CITRATED WHOLE BLOOD SPECIFICATION: 1/5/10/20/25/30/40/50 TEST/KIT

# Automated PCR analyzing system

**04**Molecular diagnostics

Real-time PCR system /14
Feature Molecular Detection Kits/14



· Available for DNA extraction and RT-PCR

 $\cdot$  Capable for testing, analyzing and screening, PC free



- ·Suitable for most of PCR systems and agents
- ·Stable in continuous experiments

Applicable reaction system	15 – 1500 μL	
Available reaction tunnels	96	
Repeatability	CV≤3.00%	
linearity	r ≥0.99	
Detecting wavelength	(500-800) nm	

# Feature Molecular Detection Kits

- ·2019-nCoV RT-PCR Test kit
- ·2019-nCoV /FluA/FluB RT-PCR Test Kit
- ·2019-nCoV/FluA/FluB/RSV/ADV RT-PCR Test Kit
- · Respiratory pathogens six-fold detection kit (PCR fluorescent probe method)
- · Monkeypox virus Nucleic Acid Diagnostic Kit (PCR-Fluorescence Probing)





Pro-med (Beijing) Technology Co., Ltd